



KUNSHAN DYNMED MEDICAL TECHNOLOGY CO. , LTD

Add: Floor 2 of No. 2 Workshop, No. 66, Shunxing Road, Lujia Town, Kunshan, Jiangsu Province, China

TEL: +86-512-82612806 FAX: +86-512-82612806

www.dynmed.cn



Product Name	Medical Oxygen Concentrator
Model Name	DO2-5E
Oxygen Flow	5L/PM
Oxygen Purity	93+/-3%
Outlet Pressure	130Kpa-150Kpa
Flow meter	Single flow
Time Counter	Yes
Alarms	Power failure, high/low pressure/Oxygen low/Working fault
Working Noise	40dB
Net Weight	13.5kg
Wheels	Removable 4 wheels
Power supply	220-240V/ 50Hz
Power consumption	300W
Power cable	European/US/UK
Electrical safety	Protection class II, type B
Filters provided	Inlet filter, bacterial filter and Coarse filter